








NEWCLIP
TECHNICS



ACTIV
MOTION

HIGH TIBIAL OSTEOTOMY &
DISTAL FEMORAL OSTEOTOMY PLATES

ACTIVMOTION S

Indications: The ACTIVMOTION S range is intended for knee osteotomy in adults


Contraindications:

- Serious vascular deterioration, bone devitalization.
- Pregnancy.
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone quality preventing the correct insertion of the screws into the bone.
- Muscular deficit, neurological deficiency or behavioural disorders, which could submit the implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

OPENING WEDGE HIGH TIBIAL OSTEOTOMY PLATES SIZES 1 & 2

→ TECHNICAL FEATURES

- ▶ **Anatomic asymmetrical implants** (blue anodized for left plates and green anodized for right plates).
- ▶ **Antero-medial positioning** to be as close as possible to the lateral stress.
- ▶ The design of the size 2 implant is adapted to large biplanar cuts or large osteotomies.
- ▶ Titanium alloy TA6V implants for optimized mechanical resistance.



→ SIZE CHOICE

▶ SIZE 1


For monoplanar osteotomy

- Up to 12 mm of correction

For biplanar osteotomy

- Up to 6 mm of correction

Ref: ATDP1D




▶ SIZE 2

For biplanar osteotomy

- Ascendant osteotomy
- Over 12 mm of correction

Ref: ATDP2D




6 to 8 monoaxial locking holes (Oneclip®) for an optimized stability

PLATE FEATURES

OPENING WEDGE HIGH TIBIAL OSTEOTOMY WITH ACL REPLACEMENT PLATES

→ TECHNICAL FEATURES

- **Anatomic asymmetrical implants** (blue anodized for left plates and green anodized for right plates).
- To prevent any risk of damaging the tunnel, the plate's upper part is **optimized for ACL reconstruction (1)**.
- 1 polyaxial locking hole located in the proximal part of the ACL tunnel to avoid damaging the graft.
- One design compatible with the peek or titanium endobutton placement.



Ref: ALTDP1D

Ref: AETDP1D






PLATE FEATURES

DISTAL FEMORAL OSTEOTOMY PLATES

→ TECHNICAL FEATURES

- **Medial closing and lateral opening plates.**
- **Anatomic asymmetrical implants** (green anodized for right plates and blue anodized for left plates).
- 2 offset screw holes improving the mechanical features of the assembly and preventing loss of angular correction (a):
 - On both sides of the osteotomy site for closing;
 - Above the osteotomy site for opening.
- Monoaxial locking screws (Oneclip®):
 - 7 screws for the closing plate;
 - 8 screws for the opening plate.
- 1 polyaxial locking screw (DTS) allowing to avoid the intercondylar notch, if necessary (b). Possible angulation of the screw before locking (25° locking range) thanks to the DTS system.
- 1 ramp oblong hole allowing for a simple and controlled compression (c) for closing.




TECHNICAL FEATURES

FIXATION FEATURES

→ SCREW TECHNICAL FEATURES


- Ø4.5 mm reinforced core screws for optimized mechanical stability (1).
- Self-tapping systems to help for insertion (2).
- Hexalobular T20 (3).



→ LOCKING SYSTEM FEATURE

➤ Low profile construct:

- The screw is stopped in the hole by its cap, insuring the locking (4).
- The screw head is buried in the plate (5) to minimize the risk of soft tissue irritation.
- Coaptation of both profiles when locking (6).
- Plate and screws made from the same material: titanium alloy.



➤ Monoaxial locking fixation


Oneclip®: patented design

➤ Polyaxial locking fixation

The DTS system (patented design) allows the screw to lock into the plate while permitting an angulation of the screw.


Newclip Technics plates combine both polyaxial and locking technologies to create a fixed-angle construct particularly useful to preserve the articulation.

Possible angulation of the screw before locking (25° locking range) thanks to the DTS® System.



→ COMPRESSIVE RAMP OBLONG HOLE


The ramp oblong hole allows a simple and controlled compression by the screw/plate interface.






TECHNICAL FEATURES

INSTRUMENTATION

- Dedicated instruments to create and maintain the appropriate angular correction during osteosynthesis:
 - 7 metallic wedges (6 mm to 18 mm; increment of 2 mm);
 - Meary pliers (controlled opening thanks to the markings - 3 to 19 mm (2 mm increment));
 - Chisels;
 - Bone spreader;
 - Cutting guide (see page 12 to see how to use the instrument);
 - MIS instruments (ANC1063 and ANC1065).



→ DIFFERENT METHODS OF OPENING OSTEOTOMIES




SURGICAL TECHNIQUE

HIGH TIBIAL OSTEOTOMY APPROACH

Technique presented below is one of the surgical techniques possibilities. The choice is made according to surgeon's preferences. Use an antero-medial approach to expose the proximal tibia metaphysis.




1. The patient is positioned in a supine position on the operating table. The procedure is performed under pneumatic tourniquet and a small pillow is placed under the buttock of the operated side in order to maintain the limb in neutral position.
2. An 8 cm slightly oblique vertical incision is made along the antero-medial surface, running over the joint space down to under the tibial tuberosity.




3. A single-plane incision is made through the periosteum; then the hamstring and the medial collateral ligament (MCL) are retracted posteriorly.

The wider the angular correction is, the more the hamstring and MCL should be released distally.

CAUTION : if the release is adequate, the opening of the osteotomy and the insertion of the bone graft can be performed with no risk of tearing the lateral cortical hinge. If it is not, forcing the graft in may tear the hinge, thus seriously jeopardizing complete bone healing.
ie: pseudarthrosis.




4. A hohmann retractor (4450-R) is placed very carefully over the posterior surface of the tibial metaphysis and should remain in place as a protection during the osteotomy.




5. Clear the deepest part of the patellar tendon down to its attachment onto the tibial tuberosity, and protect it using a retractor during the osteotomy.

SURGICAL TECHNIQUE

HIGH TIBIAL OSTEOTOMY - MONOPLANAR CUT (PAGE 1 / 2)



1. To perform the osteotomy cut, insert:
 - **The first pin** from the insertion of the hamstring until reaching the lateral cortex, 15 mm below the tibial plateau ridge.
 - **The second pin** parallel to the first one with a 20 mm distance between both of the pins, to maintain the tibial slope.




2. Incise upwards toward the head of the fibula and stop the cut 10 mm before the lateral cortical area. Then, remove the pins.



3. Insert wedges of increasing sizes until finding the appropriate one (6 - 18 mm) while maintaining the lateral surface of the tibia. Once the appropriate wedge is inserted, the angular correction is maintained during osteosynthesis.


Alternatively, chisels, the meary pliers or the bone spreader can be used to increase the size of the opening (see page 6 for more information).



4. Position the plate onto the antero-medial side so that:
 - The proximal part of the plate runs parallel to the osteotomy cut, or
 - The distal part of the plate runs parallel to the tibial tuberosity.


SURGICAL TECHNIQUE

HIGH TIBIAL OSTEOTOMY - MONOPLANAR CUT (PAGE 2/2)



5a. Lock the first Ø4.0 mm guide (ANC998) in the hole under the osteotomy cut, then start drilling using a Ø4.0 mm drill (ANC211) (1). Above the osteotomy cut, insert a Ø4.0 mm guide into the central hole (2) and drill.

Alternatively, before drilling, the plate can be temporarily maintained in position with Ø2.2 mm pin (33.0222.200) inserted through the retractor of drill guide (ANC1009)(3).




5b. The screw length can be directly read on the drill at the rear of the drill guide (see image 5a) or thanks to the length gauge (ANC210).

N.B: to ease the insertion of the screws, use the countersink (ANC120-US) to widen the first cortex previously drilled.



ANC120-US



6. Remove the drill guides. Insert and lock the 2 Ø4.5 mm screws (ST4.5LxxD-ST) using the screwdriver (ANC975). Proceed similarly for the other 4 monoaxial locking holes.



ANC975





FINAL RESULT

The construct is complete when the metallic wedge is removed

SURGICAL TECHNIQUE

HIGH TIBIAL OSTEOTOMY - BIPLANAR CUT*

* A biplanar cut must be performed with an Activmotion plate size 2





The osteotomy cut is performed in two steps :

- 1. Ascending osteotomy cut:** the cut is performed by oscillating saw, alongside and below the two pins. Stop incision 5 - 10 mm from the lateral cortex area.
- 2. Transverse osteotomy cut:** perform the anterior transverse osteotomy cut behind the tibial tuberosity at a resulting angle of around 110° to ascending cut.

- 3. Insert wedges of increasing sizes until the appropriate one while maintaining the lateral surface of the tibia. Once the appropriate wedge is inserted, the angular correction is maintained during osteosynthesis.**

Position the plate onto the antero-medial side so that the distal part of the plate runs parallel to the tibial tuberosity.

Alternatively, chisels, the meary pliers or the bone spreader can be used to increase the size of the opening (see page 6 for more information on these techniques).



- 4. Lock the first Ø4.0 mm guide (ANC998) in the hole under the osteotomy cut, then start drilling using a Ø4.0 mm drill (ANC211) (1). Above the osteotomy cut, insert a Ø4.0 mm guide into the central hole (2) and drill. Alternatively, before drilling, the plate can be temporarily maintained in position with Ø2.2 mm pin (33.0222.200) inserted through the reductor of drill guide (ANC1009)(3).**

- 5. The screw length can be directly read on the drill at the rear of the drill guide (see image 4) or thanks to the length gauge (ANC210). Remove the drill guides. Insert and lock the two Ø4.5 mm screws (ST4.5LxD-ST). Proceed similarly for the remaining locking holes.**

N.B: to ease the insertion of the screws, use the countersink (ANC120-US) to widen the first cortex previously drilled.



ANC120-US


The construct is complete when the metallic wedge is removed

SURGICAL TECHNIQUE

HIGH TIBIAL OSTEOTOMY PLATE WITH ACL REPLACEMENT



1. Perform the ACL tunnel following the surgeon's surgical technique.



2. Perform the osteotomy cut, by inserting wedges of increasing sizes until finding the appropriate one (6 - 18 mm) while maintaining the lateral surface of the tibia. Once the appropriate wedge is inserted, the angular correction is maintained during osteosynthesis.


Alternatively, chisels, the meary pliers or the bone spreader can be used to increase the size of the opening (see page 6 for more information).




3. Insert the spacer (Ø8 mm: ANC649 or Ø10 mm: ANC601) in order to preserve the tunnel during the insertion of the proximal screws (see steps 5 and 6).



4. Position the plate: the diaphyseal part of the implant should run alongside the anterior tibial tuberosity, the anterior holes are positioned on either sides of the tunnel.



5. Insert the Ø4.5 screws (ST4.5LxD-ST) located on both sides of the osteotomy site. Drill with a Ø4.0 mm drill bit (ANC211) using the drill guide (ANC998). To avoid drilling through the tunnel, use the polyaxiality for the placement of the screw into the proximal central hole. Before drilling, a pin can be inserted through the reductor of the drill guide (ANC1009)(1).




6. Once the first two screws have been inserted, repeat the procedure with the other two proximal Ø4.5 mm screws.



FINAL RESULT

Complete the procedure by inserting the last two distal screws and removing the metallic wedge and the spacer. The ligamentoplasty can then be performed.



OPTIONAL

If the plate compatible with peek or titanium endobutton is used, the endobutton (L12 x 4 mm) is inserted into the dedicated counter-forme.


SURGICAL TECHNIQUE

MEDIAL CLOSING WEDGE PLATE (PAGE 1 / 2)



1. Perform the osteotomy:

Start the first cut approximately at 5 mm above the patella groove. The cut should end around 10 mm from the lateral cortical bone.



2. Perform the proximal 2nd cut:

- Set the chosen correction angle on the cutting guide (ANC014-1 / ANC014-2).
- Insert the blade of the cutting guide into the first cut.
- Perform the second osteotomy in the cutting slot with an oscillating saw.


Remove the bone wedge and make sure that every residual bone fragment has been removed from the osteotomy.

Then, carefully close the osteotomy by applying continuous pressure to the lateral lower limb while stabilizing the knee joint region.



3. Position the plate onto the medial surface of the distal femur. The polyaxial hole must be positioned around 1 cm above the insertion of the medial collateral ligament.

It is important to ensure that the bridge area of the plate is located onto the osteotomy site and that the distal screws do not penetrate the joint.



4. Lock the first Ø4.0 mm drill guide (ANC998) into the hole situated below the osteotomy cut, drill using Ø4.0 mm drill bit (ANC211). Determine the screw length directly on the drill (1), at the rear of the drill guide or with the length gauge (ANC210) (2).



Then, insert the Ø4.5 mm locking screw (ST4.5LxxD-ST) using the screwdriver (ANC975).

Repeat this procedure with the 2 other distal monoaxial holes.

N.B. To ease the insertion of the Ø4.5 mm locking screw, use the countersink (ANC120-US) to widen the first cortex previously drilled.




How to use the cutting guide:



1. Place the NCT cutting guide - piece 2 (2) into the NCT cutting guide - piece 1 (1), and screw on the handle for metallic wedge (3).



2. Choose the correct side: **R** for right and **L** for left




3. Choose the correct angle by sliding the handle in a vertical movement; once in the correct position, turn the handle to fix in place.



4. The blade can then be inserted into the top slot of the cutting guide to perform the cut.

SURGICAL TECHNIQUE

MEDIAL CLOSING WEDGE PLATE (PAGE 2/2)



5. Insert a $\varnothing 2.2$ mm pin (33.0222.200) into the **proximal part of the oblong hole for pin**. Drill into the **proximal part of the ramp oblong hole** using the dedicated drill guide (ANC1064) and the $\varnothing 3.5$ mm drill bit (ANC1075). The orientation of the drill guide must be taken into account to allow compression (1).

Determine the screw length directly on the drill at the rear of the drill guide or with the length gauge (ANC210).



Insert a $\varnothing 4.5$ mm standard cortical screw (CT4.5LxxD-ST) and perform compression using the screwdriver (ANC975).



Then remove the $\varnothing 2.2$ mm pin.

6. Proceed similarly to step 4 for the insertion of the $\varnothing 4.5$ mm locking screws (ST4.5LxxD-ST) into the 2 holes situated above the osteotomy cut.

7. Lock the $\varnothing 4.0$ mm drill guide (ANC998) into the polyaxial hole. If necessary, adjust the drilling direction in order to avoid the intercondylar notch. Before drilling, a pin can be inserted through the reductor of the drill guide (ANC1009) (1). Start drilling using the $\varnothing 4.0$ mm drill bit (ANC211). Determine the screw length directly on the drill, at the rear of the drill guide or with the length guide (ANC210).




Then, insert the $\varnothing 4.5$ mm locking screw (ST4.5LxxD-ST) using the screwdriver (ANC975).



Repeat the previous steps to insert the remaining $\varnothing 4.5$ mm locking screws situated on the proximal part of the plate.

SURGICAL TECHNIQUE

LATERAL OPENING WEDGE PLATE (PAGE 1/2)



1. Perform the osteotomy using an oscillating saw: the cut starts 30 mm from the insertion of the lateral ligament and ends at around 10 mm from medial cortex.




2. Insert wedges of increasing sizes until finding the appropriate one (6 - 18 mm) while maintaining the lateral surface of the femur. Once the appropriate wedge has been inserted, the angular correction is maintained during osteosynthesis.


Alternatively, chisels or the meary pliers can be used to increase the size of the opening (see page 6 for more information on these techniques).



3. Position the plate onto the lateral surface of the distal femur. The polyaxial hole must be positioned at the level of the insertion of the lateral collateral ligament.



4. Stabilize the plate using the pins (33.0222.200). Insert the first pin on the distal part (1). Then, insert the second pin on the proximal part (2) so that the plate is placed alongside the femoral diaphysis.



5. Lock the first Ø4.0 mm drill guide (ANC998) in the hole under the osteotomy cut, then start drilling using the Ø4.0 mm drill bit (ANC211). Remove the drill guide and then insert and lock the Ø4.5 mm screw (ST4.5LxxD-ST) the screw using the screwdriver (ANC975)


N.B: to ease the insertion of the screws, use the countersink (ANC120-US) to widen the first cortex previously drilled.




ANC120-US

SURGICAL TECHNIQUE

LATERAL OPENING WEDGE PLATE (PAGE 2/2)



6. Proceed similarly for the insertion of the screw into the hole situated above the osteotomy cut.
Then, remove the two pins.



7. Lock the $\varnothing 4.0$ mm drill guide (ANC998) into the polyaxial hole. If necessary, adjust the drilling direction in order to avoid the intercondylar notch. Before drilling, a pin can be inserted through the reductor of the drill guide (ANC1009)(1). Start drilling using the $\varnothing 4.0$ mm drill bit (ANC211). Determine the screw length directly on the drill, at the rear of the drill guide or with the length gauge (ANC210).



ANC210

Then, insert the $\varnothing 4.5$ mm locking screw (ST4.5LxxD-ST) using the screwdriver (ANC975).

Remove the drill guide and then insert the screw using the screwdriver (ANC975)



FINAL RESULT

Repeat previous steps to insert the remaining $\varnothing 4.5$ mm locking screws. The construct is complete when the metallic wedge is removed.

IMPLANT REFERENCES


OPENING WEDGE TIBIAL PLATES

| Ref. | Description |
|------------|--|
| ATGP1D-ST | Medial opening wedge HTO plate - Left - Size 1 - STERILE |
| ATDP1D-ST | Medial opening wedge HTO plate - Right - Size 1 - STERILE |
| ATGP2D-ST | Medial opening wedge HTO plate - Left - Size 2 - STERILE |
| ATDP2D-ST | Medial opening wedge HTO plate - Right - Size 2 - STERILE |
| ALTGP1D-ST | Medial opening wedge HTO plate with ACL replacement - Left - STERILE |
| ALDTP1D-ST | Medial opening wedge HTO plate with ACL replacement - Right - STERILE |
| AETGP1D-ST | Medial opening wedge HTO plate with ACL replacement and endobutton - Left - STERILE |
| AETDP1D-ST | Medial opening wedge HTO plate with ACL replacement and endobutton - Right - STERILE |



DISTAL FEMORAL PLATES


| Ref. | Description |
|------------|--|
| JATGL1D-ST | Lateral opening wedge DFO plate - Left - Size 1 - STERILE |
| JATDL1D-ST | Lateral opening wedge DFO plate - Right - Size 1 - STERILE |
| JBTGM1D-ST | Medial closing wedge DFO plate - Left - Size 1 - STERILE |
| JBTDM1D-ST | Medial closing wedge DFO plate - Right - Size 1 - STERILE |



Ø4.5 MM DTS® SELF TAPPING SCREWS*

| Ref. | Description |
|--------------|--|
| ST4.5LxxD-ST | DTS® self-tapping screw Ø4.5 mm - L15 to 90 mm (3 mm increments from 15 to 30 mm) (5 mm increments from 30 to 90 mm) |

* Blue anodized



Ø4.5 MM CORTICAL SCREWS*

| Ref. | Description |
|--------------|--|
| CT4.5LxxD-ST | Standard cortical screw - Ø4.5 mm - L30 to 60 mm - STERILE (5 mm increments) |

* Not anodized

INSTRUMENT REFERENCES

| INSTRUMENTS | | |
|-------------|---|-----|
| Ref. | Description | Qty |
| ANC014-1 | NCT Cutting guide – piece 1 | 1 |
| ANC014-2 | NCT Cutting guide – piece 2 | 1 |
| ANC019 | Metallic wedge for knee osteotomy - 6 mm high | 1 |
| ANC020 | Metallic wedge for knee osteotomy - 8 mm high | 1 |
| ANC021 | Metallic wedge for knee osteotomy - 10 mm high | 1 |
| ANC022 | Metallic wedge for knee osteotomy - 12 mm high | 1 |
| ANC023 | Metallic wedge for knee osteotomy - 14 mm high | 1 |
| ANC024 | Handle for metallic wedge and cutting guide | 2 |
| ANC025 | Metallic wedge for knee osteotomy - 16 mm high | 1 |
| ANC120-US | Ø4.2 mm countersink with US quick coupling system | 1 |
| ANC210 | Length gauge for Ø4.5 mm screws | 1 |
| ANC211 | Ø4.0 mm quick coupling drill bit | 2 |
| ANC352 | Ø6 mm US quick coupling handle | 2 |
| ANC601 | Spacer Ø10 mm for HTO with ACL replacement | 1 |
| ANC621 | Chisel Pauwels - 10*240 mm | 1 |
| ANC622 | Chisel Pauwels - 25*240 mm | 1 |
| ANC623 | Chisel Pauwels - 35*240 mm | 1 |
| ANC628 | Chisel Pauwels - 15*240 mm | 1 |
| ANC629 | Chisel Pauwels - 20*240 mm | 1 |
| ANC649 | Spacer Ø8 mm for HTO with ACL replacement | 1 |
| ANC860 | Metallic wedge for knee osteotomy - 18 mm high | 1 |
| ANC975 | T20 screwdriver with US quick coupling system | 2 |
| ANC980 | T20 screwdriver with AO quick coupling system | 1 |
| ANC989 | Bone spreader | 1 |
| ANC990 | Activmotion Meary pliers | 1 |
| ANC998 | Ø4.0 mm threaded guide gauge | 2 |
| ANC1009 | Reducer of drill guide for Ø2.2 mm pin | 2 |
| ANC1063* | Ø4 mm mini invasive threaded guide | 2 |
| ANC1064 | Ø3.5 mm non threaded bent guide gauge | 1 |
| ANC1065* | Ø4.0 mm quick coupling drill bit - L 225 mm | 2 |
| ANC1066 | Activmotion Meary pliers | 1 |
| ANC1075 | Ø3.5 mm quick coupling drill bit - L195 mm | 2 |
| 33.0222.200 | Pin Ø2.2 L 200 mm | 6 |
| 4550-R | Hohmann retractor radiolucent | 1 |

* These instruments are dedicated to MIS

| OPTIONAL INSTRUMENTS | | |
|----------------------|-------------------------------|-----|
| Ref. | Description | Qty |
| ANC652 | HTO Alignment rod | 3 |
| ANC653 | Support for HTO alignment rod | 1 |

Tibial osteotomy plates are anatomic implants. **Bending pliers** are available on request for complex knee osteotomy cases (ANC240). Please refer to the IFU for bending precautions.

| BONE SUBSTITUTES | | |
|------------------|---------------------|-----|
| Ref. | Description | Qty |
| 0106C01 | Rounded wedge 06 mm | 1 |
| 0108C01 | Rounded wedge 08 mm | 1 |
| 0110C01 | Rounded wedge 10 mm | 1 |
| 0112C01 | Rounded wedge 12 mm | 1 |
| 1414C01 | Rounded wedge 14 mm | 1 |

Manufacturer : BIOMATLANTE (FRANCE)
Class : III
Notified body: TUV - CE 0123





REMOVAL KIT

If you have to remove ACTIVMOTION S implants, make sure to order the **Newclip Technics removal set** which includes the following instruments:
- ANC975 : T20 screwdriver with US quick coupling system
- ANC352: Ø6 mm US quick coupling handle

An **extraction set** can also be ordered separately.


Patient Specific Instruments (PSI) are also available. For more information, please refer to the **Activmotion-PSI** brochure

KIT DESCRIPTION




CLINICAL CASES

→ CASE 1: OPENING WEDGE SIZE 2




Preoperative xray




Postoperative xray

→ CASE 2: HTO & ACL



Preoperative xray



Postoperative xray

The information presented in this brochure is intended to demonstrate a NEWCLIP TECHNICS product. Always refer to the package insert, product label and/or user instructions before using any NEWCLIP TECHNICS product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your NEWCLIP TECHNICS representative if you have questions about the availability of NEWCLIP TECHNICS products in your area.



NEWCLIP TECHNICS

.....
INNOVATION MEANS MOTION

NEWCLIP TECHNICS

PA de la Lande Saint Martin
45 rue des Garottières
44115 Haute Goulaine, France
Phone: +33 (0)2 28 21 37 12
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP TECHNICS USA

Newclip USA
642 Larkfield Center
Santa Rosa CA 95403, USA
Phone: +1 707 230 5078
customerservice@newclipusa.com
www.newclipusa.com

NEWCLIP TECHNICS GERMANY

Newclip GmbH
Pröllstraße 11
86157 Augsburg, Deutschland
Phone: +49 (0)821 650 749 40
info@newclipgmbh.com
www.newclipgmbh.com

NEWCLIP TECHNICS JAPAN

Newclip Technics Japan K.K.
KKK Bldg. 502, 3-18-1 Asakusabashi
Taito-Ku, Tokyo, 111-0053, Japan
Phone: +81 z(0)3 58 25 49 81
Fax: +81 (0)3 58 25 49 86
www.newcliptechnics.com

NEWCLIP TECHNICS AUSTRALIA

Newclip Australia
3B/11 Donkin Street
West End 4101, Australia
Phone: +61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com