


INNOVATION MEANS MOTION



ALIANS PROXIMAL HUMERUS

POLYAXIAL LOCKING SYSTEM

DualTec System ®

- Polyaxiality of 25°
- Unique suture holes for soft tissues fixation
- Precontoured implant

ALIANS PROXIMAL HUMERUS

Indications: The Alians Proximal Humerus range is intended for the fixation of fractures and fracture-dislocations, osteotomies, and non-unions of the proximal humerus.

Contra-indications:

- Serious vascular deterioration, bone devitalization.
- Pregnancy.
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency touching the focus.
- Insufficient bone quality preventing the correct insertion of the screws into the bone.
- Muscular deficit, neurological deficiency or behavioral disorders, which could submit the implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- · Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

TECHNICAL FEATURES

ANATOMICALLY SHAPED PLATE

→ IDEAL PLACEMENT OF THE PLATE

- 1.5 cm from the proximal edge of the greater tuberosity, avoiding any problem of impingement.
- Alongside the bicipital groove.




- Unique & easily accessible suture holes allowing for effective stabilization of the tuberosities.
- Optimized oblong hole length allowing for adjustment of plate height.

SCREW FIXATION FEATURES

→ BLUNT-TIPPED SCREWS

- Limit protrusion through the articular surface.
- Allow to be as close as possible to the articular surface for a better construct.




Screw diameter: 4.5 mm Core diameter: 3.5 mm

TECHNICAL FEATURES

→ DEDICATED CENTERING SCREW HOLE

The drill guide (ANC131) and reductor (ANC147), with the Ø2.0 mm pin (33.0220.210), ensure an accurate placement of a centering screw in the humeral head (hole #1).

The centering screw determines the best plate positioning and optimizes the fixed-angle screw placement.





→ OPTIMIZED PROXIMAL SCREW POSITION


Divergent fixed-angle screws (targeting the inferior quadrants) and polyaxial locking screws allowing for optimized position in humeral head.





→ POLYAXIAL LOCKING SCREWS

Three variable angle locking screws allow to cope with every fracture pattern.





SURGICAL TECHNIQUE

STEP 1

SURGICAL APPROACH












The patient is placed in the beach-chair position.

A deltopectoral approach, passing outside of the cephalic vein, is recommended.

Retract the cephalic vein laterally and the pectoralis major medially.

STEP 2

FRACTURE REDUCTION



Reduce the fracture through traction and manipulation and provisionally stabilize the fracture fragments with pins (33.0220.210).

In valgus fracture patterns, the head must be elevated prior to provisional fixation.

The greater tuberosity is anatomically reduced and pinned to the shaft.


This is facilitated by manipulating the tuberosity with sutures placed through the substance of the infraspinatus. These sutures will later be used as supplemental fixation when they are secured to the plate.

Image intensification is necessary to confirm reduction.

STEP 3

OSTEOSYNTHESIS PROCEDURE

→ CENTERING SCREW



Place the plate alongside the bicipital groove and approximately 1.5 cm distal to the top of the greater tuberosity.

Insert a \emptyset 4.5 mm cortical screw (CT4.5Lxx) into the oblong hole and fasten the plate to the shaft. Provisionally secure the plate to the bone with \emptyset 2.0 mm pins (33.0220.210).

Insert the drill guide (ANC131) with its reductor (ANC147) through hole #1. Insert a \emptyset 2.0 mm pin (33.0220.210) to target the center of the humeral head. Check position and trajectory with the C-arm.

Then drill at \emptyset 3.5 mm (ANC132) through the drill guide (ANC131) and insert the first \emptyset 4.5 mm locking screw (PT4.5Lxx).



SURGICAL TECHNIQUE

→ FIXED-ANGLE DIVERGENT SCREWS




Use the $\emptyset 3.5$ mm drill guide (ANC131) and drill 4 divergent $\emptyset 4.5$ mm fixed-angle screws (PT4.5Lxx).

This precise screw pattern enhances resistance to varus forces.

Blunt-tipped screws limit protrusion through the articular surface.

→ POLYAXIAL LOCKING SCREWS



Orientate and lock the first 2 proximal screws and the first metaphyseal screw (1) according to the fracture pattern.

As the highest bone density is located in the inferior quadrants, every attempt should be made to keep the screws descending.

Use the $\emptyset 3.5$ mm drill guide (ANC127) and place the remaining distal cortical screws, non locking (CT4.5Lxx) or locking (VT4.5Lxx) at the surgeon's preference.



STEP 4 SUTURE OF THE TUBEROSITIES AND C-ARM CONTROL




Repair and fasten the tuberosity to the plate through the suture holes.


Assess the final reduction under fluoroscopy.

IMPLANTS REFERENCES

- RIGHT PLATES - GREEN ANODIZED




→ LEFT PLATES - BLUE ANODIZED



	ALIANS PROXIMAL HUMERUS PLATES
Ref.	Description
STDPSS1	Proximal humerus plate - Right - Size 1 - Short
STGPSS1	Proximal humerus plate - Left - Size 1 - Short
STDPS1	Proximal humerus plate - standard - Right - Size 1
STGPS1	Proximal humerus plate - standard - Left - Size 1
STDPS2	Proximal humerus plate - standard - Right - Size 2
STGPS2	Proximal humerus plate - standard - Left - Size 2
STDPS3	Proximal humerus plate - standard - Right - Size 3
STGPS3	Proximal humerus plate - standard - Left - Size 3

Suture holes for soft tissue fixation are compatible with Ø2.0 mm needles.

IMPLANTS REFERENCES



^{*} Green anodized

The state of the s	Ø4.5 MM
	NON LOCKED RETAINING SCREWS*
Ref.	Description
QT4.5L32	Non locked retaining screws - Ø4.5 mm - L32 mm
QT4.5L36	Non locked retaining screws - Ø4.5 mm - L36 mm
QT4.5L40	Non locked retaining screws - Ø4.5 mm - L40 mm
QT4.5L44	Non locked retaining screws - Ø4.5 mm - L44 mm

^{*} Golden anodized

Ø4.5 MM STANDARD CORTICAL SCREWS* Standard cortical screw - Ø4.5 mm - L20 mm CT4.5L20 CT4.5L22 Standard cortical screw - Ø4.5 mm - L22 mm CT4.5L24 Standard cortical screw - Ø4.5 mm - L24 mm Standard cortical screw - Ø4.5 mm - L26 mm CT4.5L26 Standard cortical screw - Ø4.5 mm - L28 mm CT4.5L28 CT4.5L30 Standard cortical screw - Ø4.5 mm - L30 mm CT4.5L32 Standard cortical screw - Ø4.5 mm - L32 mm CT4.5L34 Standard cortical screw - Ø4.5 mm - L34 mm CT4.5L36 Standard cortical screw - Ø4.5 mm - L36 mm CT4.5L38 Standard cortical screw - Ø4.5 mm - L38 mm CT4.5L40 Standard cortical screw - Ø4.5 mm - L40 mm

* Not anodized

4	
N.	Ø4.5 MM
~	LOCKING SELF-TAPPING CORTICAL SCREWS*
Ref.	Description
VT4.5L20	Locking self-tapping cortical screw - Ø4.5 mm - L20 mm
VT4.5L22	Locking self-tapping cortical screw - Ø4.5 mm - L22 mm
VT4.5L24	Locking self-tapping cortical screw - Ø4.5 mm - L24 mm
VT4.5L26	Locking self-tapping cortical screw - Ø4.5 mm - L26 mm
VT4.5L28	Locking self-tapping cortical screw - Ø4.5 mm - L28 mm
VT4.5L30	Locking self-tapping cortical screw - Ø4.5 mm - L30 mm
VT4.5L32	Locking self-tapping cortical screw - Ø4.5 mm - L32 mm
VT4.5L34	Locking self-tapping cortical screw - Ø4.5 mm - L34 mm
VT4.5L36	Locking self-tapping cortical screw - Ø4.5 mm - L36 mm
VT4.5L38	Locking self-tapping cortical screw - Ø4.5 mm - L38 mm
VT4.5L40 * Blue anodized	Locking self-tapping cortical screw - Ø4.5 mm - L40 mm

Remark:



Please note that all implants are also available in sterile packaging. The double tube packaging is handy and easy to use. An "ST" code is added at the end of the reference. e.g. "PT4.5L26-ST"

RACK (ANC133/R)



(ANC133/B)

INSTRUMENTS				
Ref.	Description	Qty		
ANC119-US	3.0 mm quick coupling hexagonal screwdriver	2		
ANC120-US	Ø4.2 mm reamer with US quick coupling system	1		
ANC121-US	Ø3.5 mm drill guide with US quick coupling system	1		
ANC127	Ø3.5 mm drill guide for locking cortical screws Ø4.5 mm	2		
ANC129	Length gauge for Ø4.5 mm DTS® screws	1		
ANC131	Ø3.5 mm drill guide for DTS® screws	2		
ANC132	Ø3.5 mm quick coupling drill bit - L195 mm	2		
ANC147	Reductor of drill guide DTS Shoulder for Ø2.0 mm pin	1		
ANC352	Ø6 mm US quick coupling handle	2		
33.0220.210	Pin Ø2.0 L210 mm	3		

REMOVAL SET

If you have to remove Alians Proximal Humerus implants, make sure to order the Newclip Technics removal set which includes the following

- ANC119-US: 3.0 mm quick coupling hexagonal screwdriver for Ø4.5 mm screws
- ANC352: Ø6 mm US quick coupling handle

The information presented in this brochure is intended to demonstrate a NEWCLIP TECHNICS product. Always refer to the package insert, product label and/or user instructions before using any NEWCLIP TECHNICS product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your NEWCLIP TECHNICS representative if you have questions about the availability of NEWCLIP TECHNICS products in your area.

NEWCLIP USA

