


TOUCH®

Dual mobility
Trapeziometacarpal
prosthesis

Highly Cross-linked Vs. Standard PE*



Test* on 5 million cycles
(approx. 10 years of regular thumb usage)

Highly Cross-linked Liner




*Internal testing at KERI MEDICAL SA (110-30A102)
Note : bench testing is not necessarily indicative of clinical performance

Fully committed to **Hand & Wrist** only
www.kerimedical.com

/ SPHERICAL CUPS

- 2 cup sizes (Ø 9 mm and 10 mm)
- 5 anti-rotation fins
- Spiked crown for stability
- Material: Stainless steel
- Dual coating: Titanium + HA



CTO09

CTO10

/ CONICAL CUPS

- 2 cup sizes (Ø 9 mm and 10 mm)
- Material: Stainless steel
- Dual coating: Titanium + HA



CTO109

CTO110

/ LINER AND STRAIGHT NECKS

- Liner material: Highly Cross-linked PE
- Neck material: Stainless Steel
- Straight neck
- 3 heights: 6, 8 and 10 mm
- Head diameter 4 mm



NTO06

NTO08

NTO010

/ LINER AND OFFSET NECKS

- Liner material: Highly Cross-linked PE
- Neck material: Stainless Steel
- 15° offset neck
- 3 heights: 6, 8 and 10 mm
- Head diameter 4 mm



NTO156

NTO158

NTO1510

/ STEMS

- 6 stems
- Anatomical trilobe section
- Material: Titanium
- Dual coating: Titanium + HA



STOXS

STO0

STO1

STO2

STO3

STO4

The + TOUCH®

1st

Dual mobility prosthesis with highly Cross-linked Liner


1st

prosthesis on the market to propose a spherical and conical cup

144

possible combinations to fit as close to the patient's anatomy

/ INSTRUMENT SET



+ TOUCH®

Complete set including orthopaedic standard instruments

CE 2862

KERI MEDICAL SA

Chemin du Pré-Fleuri 5
1228 Plan-les-Ouates - Switzerland
Tel: +41 58 255 01 30
Fax: +41 58 255 01 40
customerservice@kerimedical.com

EC REP KERI MEDICAL FRANCE

34 rue Antoine Redier
74160 Archamps - France
Tel: +33(0)4 50 85 36 10
Fax: +33 (0)4 80 16 05 10
adv@kerimedical.com

Distributed by :